

Product datasheet for **TA356013**

Zbtb43 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat Homology: Cow: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 93%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	57kDa
Gene Name:	zinc finger and BTB domain containing 43
Database Link:	NP_082223 Entrez Gene 71834 Mouse Q9DAI4
Background:	Zbtb43 may be involved in transcriptional regulation.
Synonyms:	FLJ22470; KIAA0414; OTTHUMP00000022157; ZBTB22B; ZNF-X; ZNF297B



[View online »](#)

Product images:



WB Suggested Anti-Zbtb43 Antibody
Titration: 1.0 ug/ml
Positive Control: Mouse Kidney