

## Product datasheet for **TA356024**

### ZNF96 (ZSCAN12) Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N-terminal region of Human ZSC12
Specificity:	<b>Expected reactivity:</b> Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	66kDa
Gene Name:	zinc finger and SCAN domain containing 12
Database Link:	<a href="#">NP_055539</a> <a href="#">Entrez Gene 9753 Human</a> <a href="#">O43309</a>
Synonyms:	dj29K1.2; KIAA0426; OTTHUMP00000016205; ZFP96; ZNF29K1; ZNF96; ZNF305



[View online »](#)

## Product images:



Host: Rabbit  
Target Name: ZSC12  
Sample Tissue: RPMI-8226 Cell Lysate  
Antibody Dilution: 1.0 $\mu$ g/ml

Host: Rabbit  
Target Name: ZSC12  
Sample Type: RPMI-8226 Whole Cell lysates  
Antibody Dilution: 1.0ug/ml