

## Product datasheet for **TA356074**

### Kat5 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Dog, Human, Mouse, Rat <b>Homology:</b> Dog: 93%; Human: 100%; Mouse: 100%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	47kDa
Gene Name:	K(lysine) acetyltransferase 5
Database Link:	<a href="#">NP_001186178</a> <a href="#">Entrez Gene 81601 Mouse</a> <a href="#">Q8CHK4</a>
Background:	The function of this protein remains unknown.
Synonyms:	cPLA2; ESA1; HTATIP; HTATIP1; PLIP; TIP; TIP60



[View online »](#)

## Product images:



WB Suggested Anti-Kat5 Antibody  
Titration: 1.0 ug/ml  
Positive Control: Mouse Thymus