

Product datasheet for **TA356085**

Rbm3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Yeast Homology: Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%; Yeast: 77%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	17kDa
Gene Name:	RNA binding motif (RNP1, RRM) protein 3
Database Link:	NP_446148 Entrez Gene 114488 Rat Q925G0-2
Background:	Mouse homolog displays increased protein synthesis in response to mild hypothermia.
Synonyms:	IS1-RNPL; OTTHUMP00000025802; RNPL



[View online »](#)

Product images:

WB Suggested Anti-Rbm3 Antibody
Titration: 1.0 ug/ml
Positive Control: Rat Lung