

Product datasheet for **TA356094**

Srp54a Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Sheep, Zebrafish Homology: Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%; Sheep: 83%; Zebrafish: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	55kDa
Gene Name:	signal recognition particle 54A
Database Link:	NP_446323 Entrez Gene 116650 Rat Q6AYB5
Background:	may recognize signal sequences of secretory proteins and may retard them in the endoplasmic reticulum; may regulate synthesis of the secretory protein osteopontin [RGD, Feb 2006]
Synonyms:	54kDa; Srp54



[View online »](#)

Product images:



WB Suggested Anti-Srp54a Antibody
Titration: 1.0 ug/ml
Positive Control: Rat Muscle