

Product datasheet for **TA356141**

Hao2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Horse, Human, Mouse, Pig, Rat Homology: Cow: 93%; Dog: 93%; Horse: 100%; Human: 100%; Mouse: 93%; Pig: 93%; Rat: 93%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	39kDa
Gene Name:	hydroxyacid oxidase 2
Database Link:	NP_062418 Entrez Gene 56185 Mouse Q9NYQ2
Background:	Hao2 has 2-hydroxyacid oxidase activity and most active on medium-chain substrates.
Synonyms:	GIG16; HAOX2



[View online »](#)

Product images:



WB Suggested Anti-Hao2 Antibody
Titration: 1.0 ug/ml
Positive Control: Mouse Brain