

## Product datasheet for **TA356146**

### Adamts4 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Yeast <b>Homology:</b> Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%; Yeast: 79%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	93kDa
Gene Name:	ADAM metalloproteinase with thrombospondin type 1 motif, 4
Database Link:	<a href="#">NP_076449</a> <a href="#">Entrez Gene 66015 Rat</a> <a href="#">Q9ESP7</a>
Background:	Adamts4 is a disintegrin and metalloproteinase; Adamts4 may play a role in the progression of Alzheimers disease.
Synonyms:	ADAM-TS4; ADAMTS-2; ADAMTS-4; ADMP-1; aggrecanase-1; KIAA0688



[View online »](#)

Product images:



Host: Rabbit  
Target Name: Adamts4  
Antibody Dilution: 1.0ug/ml  
Sample Type: Rat Kidney