

Product datasheet for **TA356165**

Sst Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Goat, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Sheep Homology: Cow: 100%; Dog: 100%; Goat: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%; Sheep: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	13kDa
Gene Name:	somatostatin
Database Link:	NP_033241 Entrez Gene 24797 Rat P60041
Background:	Somatostatin inhibits the release of somatotropin.
Synonyms:	SMST; somatostatin; somatostatin-14; somatostatin-28



[View online »](#)

Product images:



Rabbit Anti-Sst Antibody
Catalog Number: ARP41905
Lot Number: QC12533
Lane: Rat Liver Lysate

Antibody Titration: 1.0 μ g/ml
Gel Concentration: 10-20%

Host: Rabbit
Target Name: SST
Sample Tissue: Rat Liver
Antibody Dilution: 1 μ g/ml