

## Product datasheet for **TA356175**

### Btrc Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat <b>Homology:</b> Cow: 100%; Dog: 100%; Guinea Pig: 93%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	68kDa
Gene Name:	beta-transducin repeat containing E3 ubiquitin protein ligase
Database Link:	<a href="#">NP_001007149</a> <a href="#">Entrez Gene 361765 Rat</a>
Background:	The function of this protein remains unknown.
Synonyms:	BETA-TRCP; beta-TrCP1; betaTrCP; bTrCP; bTrCP1; E3RSIkappaB; FBW1A; FBXW1; FBXW1A; FWD1; MGC4643



[View online »](#)

**Product images:**

Rabbit Anti-Btrc Antibody  
Catalog Number: ARP42159  
Lot Number: QC12927  
Lane: Rat Lung Lysate

Antibody Titration: 1.0 $\mu$ g/ml  
Gel Concentration: 12%

Host: Rabbit  
Target Name: BTRC  
Sample Tissue: Rat Lung  
Antibody Dilution: 1 $\mu$ g/ml