

Product datasheet for **TA356196**

Eva1b Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Dog, Goat, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Yeast Homology: Dog: 89%; Goat: 80%; Guinea Pig: 87%; Horse: 100%; Human: 100%; Mouse: 82%; Rabbit: 86%; Rat: 93%; Yeast: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	24kDa
Gene Name:	eva-1 homolog B (C. elegans)
Database Link:	NP_001102149 Entrez Gene 362597 Rat B2RZB3
Background:	The function of this protein remains unknown.
Synonyms:	2610027C15Rik; Fam176b; RGD1308876



[View online »](#)

Product images:



WB Suggested Anti-Fam176b Antibody
Titration: 1.0 ug/ml
Positive Control: Rat Lung