

## Product datasheet for **TA356325**

### Ubox5 Rabbit Polyclonal Antibody

#### Product data:

|                         |                                                                                                                                                                                                                            |
|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:           | Primary Antibodies                                                                                                                                                                                                         |
| Applications:           | WB                                                                                                                                                                                                                         |
| Reactivity:             | Mouse                                                                                                                                                                                                                      |
| Host:                   | Rabbit                                                                                                                                                                                                                     |
| Clonality:              | Polyclonal                                                                                                                                                                                                                 |
| Specificity:            | <b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Pig, Rabbit, Rat<br><b>Homology:</b> Cow: 100%; Dog: 100%; Guinea Pig: 93%; Horse: 100%; Human: 100%; Mouse: 93%; Pig: 100%; Rabbit: 100%; Rat: 93% |
| Formulation:            | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.<br><i>Note that this product is shipped as lyophilized powder to China customers.</i>                                    |
| Concentration:          | lot specific                                                                                                                                                                                                               |
| Purification:           | Affinity Purified                                                                                                                                                                                                          |
| Conjugation:            | Unconjugated                                                                                                                                                                                                               |
| Storage:                | For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.                                                                                    |
| Stability:              | Shelf life: one year from despatch.                                                                                                                                                                                        |
| Predicted Protein Size: | 59kDa                                                                                                                                                                                                                      |
| Gene Name:              | U box domain containing 5                                                                                                                                                                                                  |
| Database Link:          | <a href="#">NP_542129</a><br><a href="#">Entrez Gene 140629 Mouse</a><br><a href="#">Q925F4</a>                                                                                                                            |
| Background:             | The function of this protein remains unknown.                                                                                                                                                                              |
| Synonyms:               | KIAA0860; RNF37; Ubce7ip5; UIP5                                                                                                                                                                                            |



[View online »](#)

Product images:



WB Suggested Anti-Ubox5 Antibody  
Titration: 1.0 ug/ml  
Positive Control: Mouse Heart