

Product datasheet for **TA356381**

SLC45A1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	96kDa
Gene Name:	solute carrier family 45 member 1
Database Link:	XP_942788 Entrez Gene 50651 Human Q9Y2W3
Background:	The function of this protein remains unknown.
Synonyms:	DNb-5; DNB5; FLJ26758; KIAA0458; PAST-A
Protein Families:	Transmembrane



[View online »](#)

Product images:



WB Suggested Anti-SLC45A1 Antibody
Titration: 1.0 ug/ml
Positive Control: COLO205 Whole Cell