

## Product datasheet for **TA356487**

### LRCH4 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human LRCH4
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Pig, Rabbit, Rat, Zebrafish <b>Homology:</b> Cow: 93%; Dog: 100%; Guinea Pig: 93%; Horse: 100%; Human: 100%; Mouse: 100%; Pig: 100%; Rabbit: 100%; Rat: 100%; Zebrafish: 83%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	73kDa
Gene Name:	leucine-rich repeats and calponin homology (CH) domain containing 4
Database Link:	<a href="#">NP_002310</a> <a href="#">Entrez Gene 4034 Human</a> <a href="#">O75427</a>
Background:	The specific function of LRCH4 is not yet known. This gene encodes a protein that contains leucine-rich repeats (LRR) at its amino terminus and that is known to be involved in ligand binding. The carboxyl terminus may act as a membrane anchor. Identified structural elements suggest that the encoded protein resembles a receptor.
Synonyms:	FLJ40101; FLJ46315; LRN; LRRN1; LRRN4; PP14183
Protein Families:	Druggable Genome, Transmembrane



[View online >](#)

Product images:



WB Suggested Anti-LRCH4 Antibody Titration: 0.2-1 ug/ml  
ELISA Titer: 1:62500  
Positive Control: Hela cell lysate