

Product datasheet for **TA356581**

Olig1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Human, Mouse, Rabbit, Rat Homology: Cow: 100%; Dog: 100%; Guinea Pig: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	28kDa
Gene Name:	oligodendrocyte transcription factor 1
Database Link:	NP_068538 Entrez Gene 60394 Rat Q9WUQ3
Background:	This gene is a member of the oligodendrocyte lineage gene family which encode basic helix-loop-helix transcription factors. Studies in mice show that it is coexpressed with Olig2 in neural progenitors, and has a role in the development and maturation of oligodendrocytes.
Synonyms:	BHLHB6; BHLHE21; Oligo1



[View online »](#)

Product images:



WB Suggested Anti-Olig1 Antibody
Titration: 1.0 ug/ml
Positive Control: Rat Lung