

Product datasheet for **TA356586**

Sik1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Dog, Guinea Pig, Horse, Human, Mouse, Rat, Zebrafish Homology: Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rat: 100%; Zebrafish: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	85kDa
Gene Name:	salt-inducible kinase 1
Database Link:	NP_067725 Entrez Gene 59329 Rat Q9R1U5
Background:	Sik1 inhibits CREB binding and transcription of cAMP-responsive element (CRE)-dependent genes Cyp11A and Star.
Synonyms:	msk; SIK; SIK-1; SNF1LK



[View online »](#)

Product images:



WB Suggested Anti-Sik1 Antibody
Titration: 1.0 ug/ml
Positive Control: Rat Lung