

## Product datasheet for **TA356618**

### SYPL1 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Dog, Horse, Human, Pig <b>Homology:</b> Dog: 85%; Horse: 77%; Human: 100%; Pig: 77%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	28kDa
Gene Name:	synaptophysin like 1
Database Link:	<a href="#">NP_006745</a> <a href="#">Entrez Gene 6856 Human</a> <a href="#">Q16563</a>
Background:	The function of this protein remains unknown.
Synonyms:	H-SP1; pantophysin; SYPL
Protein Families:	Transmembrane



[View online »](#)

## Product images:



WB Suggested Anti-SYPL1 Antibody  
Titration: 1.0 ug/ml  
Positive Control: Fetal Thymus