

Product datasheet for TA356623

Golim4 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Zebrafish Homology: Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%; Zebrafish: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2–8°C up to 1 week. For long term storage, store at –20°C in small aliquots to prevent freeze–thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	75kDa
Gene Name:	golgi integral membrane protein 4
Database Link:	NP_001178496 Entrez Gene 310526 Rat Q5BJK8
Background:	Golim4 plays a role in endosome to Golgi protein trafficking; mediates protein transport along the late endosome–bypass pathway from the early endosome to the Golgi.
Synonyms:	GIMPC; GOLPH4; GPP130; P138



Product images:



WB Suggested Anti-Golim4 Antibody
 Titration: 1.0 ug/ml
 Positive Control: Rat Muscle