

Product datasheet for **TA356649**

MCTP2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of Human MCTP2
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Pig, Rabbit, Rat Homology: Cow: 93%; Dog: 100%; Guinea Pig: 86%; Horse: 86%; Human: 100%; Mouse: 100%; Pig: 93%; Rabbit: 93%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	51kDa
Gene Name:	multiple C2 domains, transmembrane 2
Database Link:	Entrez Gene 55784 Human Q6DN12-3
Background:	The function of this protein remains unknown.
Synonyms:	FLJ11175; FLJ33303; MGC111402



[View online »](#)

Product images:



Host: Rabbit
 Target Name: MCTP2
 Sample Type: HeLa Whole Cell lysates
 Antibody Dilution: 1.0ug/ml
 MCTP2 is supported by BioGPS gene expression data to be expressed in HeLa