

Product datasheet for **TA356670**

Xpo5 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Pig, Rabbit, Rat Homology: Cow: 100%; Dog: 92%; Guinea Pig: 92%; Horse: 92%; Human: 100%; Mouse: 100%; Pig: 100%; Rabbit: 90%; Rat: 93%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	132kDa
Gene Name:	exportin 5
Database Link:	NP_001102259 Entrez Gene 363194 Rat
Background:	The function of this protein remains unknown.
Synonyms:	Exp5; FLJ14239; FLJ32057; FLJ45606; KIAA1291; RANBP21



[View online »](#)

Product images:



WB Suggested Anti-Xpo5 Antibody
Titration: 1.0 ug/ml
Positive Control: Rat Liver