

## Product datasheet for **TA356672**

### Irx6 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat <b>Homology:</b> Cow: 92%; Dog: 92%; Guinea Pig: 92%; Horse: 85%; Human: 100%; Mouse: 86%; Rabbit: 92%; Rat: 92%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	47kDa
Gene Name:	iroquois homeobox 6
Database Link:	<a href="#">NP_001095883</a> <a href="#">Entrez Gene 307715 Rat</a> <a href="#">D4A7R4</a>
Background:	The function of this protein remains unknown.
Synonyms:	IRX-3; IRX7; IRXB3



[View online »](#)

## Product images:



WB Suggested Anti-Irx6 Antibody  
Titration: 1.0 ug/ml  
Positive Control: Rat Muscle