

## Product datasheet for TA356683

### Ankyrin repeat domain containing protein 65 (ANKRD65) Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of human hCG_20426
Specificity:	<b>Expected reactivity:</b> Dog, Human <b>Homology:</b> Dog: 92%; Human: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	42kDa
Gene Name:	ankyrin repeat domain 65
Database Link:	<a href="#">XP_497648</a> <a href="#">Entrez Gene 441869 Human</a> <a href="#">E5RJM6</a>
Background:	The function remain unknown.
Synonyms:	ANKRD65



[View online »](#)

**Product images:**



WB Suggested Anti-hCG\_20426 Antibody  
Titration: 0.2-1 ug/ml  
Positive Control: Human Muscle