

Product datasheet for **TA356687**

Antp Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Drosophila
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide corresponding to a region of Fruit fly
Specificity:	Homology: Fruit fly:100%;
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	33kDa
Gene Name:	Antennapedia
Database Link:	NP_996176 P02833
Background:	Antp is a sequence-specific transcription factor which is part of a developmental regulatory system that regulates segmental identity in the mesothorax. It provides cells with specific positional identities on the anterior-posterior axis.



[View online »](#)

Product images:



WB Suggested Anti-Antp Antibody Titration: 0.2-1 ug/ml
Positive Control: Drosophila