

## Product datasheet for **TA356704**

### **E130311K13Rik Rabbit Polyclonal Antibody**


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Pig, Rabbit, Rat <b>Homology:</b> Cow: 93%; Dog: 93%; Guinea Pig: 79%; Horse: 93%; Human: 100%; Mouse: 86%; Pig: 100%; Rabbit: 93%; Rat: 86%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	34kDa
Gene Name:	RIKEN cDNA E130311K13 gene
Database Link:	<a href="#">XP_002725693</a> <a href="#">Entrez Gene 329659 Mouse</a> <a href="#">Q8BN57</a>
Background:	The function of this protein remains unknown.
Synonyms:	E130311K13Rik



[View online »](#)

Product images:



WB Suggested Anti-E130311K13Rik Antibody  
Titration: 1.0 ug/ml  
Positive Control: Mouse Heart