

Product datasheet for **TA356766**

Sult1c2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat Homology: Cow: 86%; Dog: 100%; Guinea Pig: 93%; Horse: 100%; Human: 100%; Mouse: 86%; Rabbit: 93%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	32kDa
Gene Name:	sulfotransferase family 1C member 2
Database Link:	NP_598231 Entrez Gene 171072 Rat Q9WUW8
Background:	Sult1c2 catalyzes the sulfation of p-nitrophenol; localized to lysosomes.
Synonyms:	humSULTC2; ST1C1; ST1C2; SULT1C1



[View online »](#)

Product images:



WB Suggested Anti-Sult1c2 Antibody
Titration: 1.0 ug/ml
Positive Control: Rat Liver