

## Product datasheet for **TA356772**

### Mgat5b Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rat, Zebrafish <b>Homology:</b> Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rat: 100%; Zebrafish: 79%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	67kDa
Gene Name:	mannosyl (alpha-1,6-)-glycoprotein beta-1,6-N-acetyl-glucosaminyltransferase, isozyme B
Database Link:	<a href="#">NP_001100538</a> <a href="#">Entrez Gene 303693 Rat</a> <a href="#">D3ZRS9</a>
Background:	The function of this protein remains unknown.
Synonyms:	FLJ25132; GnT-IX; GnT-VB; hGnTVb; KIAA2008



[View online »](#)

Product images:



WB Suggested Anti-Mgat5b Antibody  
Titration: 1.0 ug/ml  
Positive Control: Rat Muscle