

## Product datasheet for **TA356776**

### Eef1akmt2 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Dog, Guinea Pig, Human, Mouse, Rabbit, Rat <b>Homology:</b> Dog: 100%; Guinea Pig: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	27kDa
Gene Name:	methyltransferase like 10
Database Link:	<a href="#">NP_082371</a> <a href="#">Entrez Gene 72096 Mouse</a> <a href="#">Q9D853</a>
Synonyms:	METTL10



[View online »](#)

Product images:



WB Suggested Anti-Mettl10 Antibody  
Titration: 1.0 ug/ml  
Positive Control: Mouse Spleen