

Product datasheet for **TA356790**

GLYATL1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of Human GLYL1
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	36kDa
Gene Name:	glycine-N-acyltransferase like 1
Database Link:	NP_542392 Entrez Gene 92292 Human Q969I3-2
Synonyms:	FLJ26507; FLJ34646; GATF-C; GNAT; MGC15397; MGC15937



[View online »](#)

Product images:



Host: Rabbit
Target Name: GLYL1
Sample Tissue: COLO205 Cell Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: GLYL1
Sample Type: COLO205 Whole Cell lysates
Antibody Dilution: 1.0ug/ml