

## Product datasheet for **TA356804**

### Tm9sf3 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N-terminal region of Rat Tm9sf3
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Zebrafish <b>Homology:</b> Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%; Zebrafish: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	67kDa
Gene Name:	transmembrane 9 superfamily member 3
Database Link:	<a href="#">XP_220013</a> <a href="#">Entrez Gene 309475 Rat</a> <a href="#">D3ZUD8</a>
Background:	The function of this protein remains unknown.
Synonyms:	EP70-P-iso; RP11-34E5.1; SMBP



[View online »](#)

Product images:



Host: Rabbit  
Target Name: Tm9sf3  
Sample Type: Rat Stomach lysates  
Antibody Dilution: 1.0ug/ml