

Product datasheet for **TA356836**

Gal3st4 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat Homology: Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	52kDa
Gene Name:	galactose-3-O-sulfotransferase 4
Database Link:	NP_001102534 Entrez Gene 498166 Rat D3ZHZ0
Background:	The function of this protein remains unknown.
Synonyms:	FLJ12116; GAL3ST-4; OTTHUMP00000206900



[View online »](#)

Product images:



WB Suggested Anti-Gal3st4 Antibody
Titration: 1.0 ug/ml
Positive Control: Rat Liver