

## Product datasheet for **TA356866**

### Tmem132d Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Zebrafish <b>Homology:</b> Cow: 93%; Dog: 93%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 93%; Rabbit: 100%; Rat: 100%; Zebrafish: 86%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	121kDa
Gene Name:	transmembrane protein 132D
Database Link:	<a href="#">NP_766473</a> <a href="#">Entrez Gene 243274 Mouse Q76HP3</a>
Background:	Tmem132d may serve as a cell-surface marker for oligodendrocyte differentiation.
Synonyms:	HBE120; KIAA1944; MGC138770; MGC138771; MOLT



[View online »](#)

Product images:



WB Suggested Anti-Tmem132d Antibody  
Titration: 1.0 ug/ml  
Positive Control: Mouse Liver