

Product datasheet for **TA356903**

ZNF891 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N terminal region of human hCG_1646157
Specificity:	Expected reactivity: Human Homology: Human: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	63kDa
Gene Name:	similar to KRAB zinc finger protein 6D
Database Link:	XP_001129759 A8MT65
Background:	The exact function of hCG_1646157 remains unknown.



[View online »](#)

Product images:


WB Suggested Anti-hCG_1646157 Antibody
 Titration: 0.2-1 ug/ml
 Positive Control: Hela cell lysate