

Product datasheet for **TA356965**


Zfp874b Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Human, Mouse, Rat Homology: Cow: 79%; Human: 100%; Mouse: 100%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	52kDa
Gene Name:	zinc finger protein 874b
Database Link:	NP_001070259 Entrez Gene 408067 Mouse D3Z0W3
Background:	The function of this protein remains unknown.
Synonyms:	9630025I21Rik; Rscan16



[View online »](#)

Product images:


WB Suggested Anti-Zfp874b Antibody
 Titration: 1.0 ug/ml
 Positive Control: Mouse Kidney