

## Product datasheet for TA356980

### C4orf46 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Horse, Human <b>Homology:</b> Horse: 86%; Human: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	12kDa
Gene Name:	chromosome 4 open reading frame 46
Database Link:	<a href="#">NP_001008394</a> <a href="#">Entrez Gene 201725 Human</a> <a href="#">Q504U0</a>
Background:	The specific function of LOC201725 is not yet known.
Synonyms:	LOC201725



[View online »](#)

## Product images:



WB Suggested Anti-C4orf46 Antibody  
Titration: 1.0 ug/ml  
Positive Control: 293T Whole Cell