

Product datasheet for **TA356988**

Pnck Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Human, Mouse, Pig, Rabbit, Rat Homology: Cow: 100%; Dog: 100%; Human: 100%; Mouse: 100%; Pig: 100%; Rabbit: 100%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	37kDa
Gene Name:	pregnancy up-regulated nonubiquitous CaM kinase
Database Link:	NP_058971 Entrez Gene 29660 Rat O70150
Background:	Pnck contains alternatively spliced isoforms CaM kinase Ibeta3 and Ibeta1; both isoforms activate CaM kinase Ialpha activity via Thr177 phosphorylation.
Synonyms:	BSTK3; CaMK1b; CaMKI-beta; FLJ50403; FLJ50549; FLJ56451; FLJ59811; MGC45419



[View online »](#)

Product images:



Host: Rabbit
Target Name: Pnck
Antibody Dilution: 1.0ug/ml
Sample Type: Rat Stomach