

Product datasheet for **TA357066**

Glt1d1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the N-terminal region of Rat Glt1d1
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Human, Mouse, Rat Homology: Cow: 93%; Dog: 93%; Guinea Pig: 93%; Human: 100%; Mouse: 100%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	38kDa
Gene Name:	glycosyltransferase 1 domain containing 1
Database Link:	XP_222147 Entrez Gene 304445 Rat D3ZUX2
Background:	The function of this protein remains unknown.
Synonyms:	FLJ31978



[View online »](#)

Product images:



Host: Rabbit
 Target Name: Glt1d1
 Sample Type: Rat Pancreas lysates
 Antibody Dilution: 1.0ug/ml