

## Product datasheet for **TA357074**

### Fam122b Rabbit Polyclonal Antibody


#### Product data:

|                         |   |
|-------------------------|---|
| Product Type:           | Primary Antibodies  |
| Applications:           | WB  |
| Reactivity:             | Mouse   |
| Host:                   | Rabbit  |
| Clonality:              | Polyclonal  |
| Specificity:            | <b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat<br><b>Homology:</b> Cow: 93%; Dog: 93%; Guinea Pig: 93%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 93%; Rat: 100% |
| Formulation:            | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.   |
| Concentration:          | lot specific  |
| Purification:           | Affinity Purified   |
| Conjugation:            | Unconjugated  |
| Storage:                | For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.   |
| Stability:              | Shelf life: one year from despatch.   |
| Predicted Protein Size: | 28kDa   |
| Gene Name:              | family with sequence similarity 122, member B   |
| Database Link:          | <a href="#">NP_084443</a><br><a href="#">Entrez Gene 78755 Mouse</a><br><a href="#">Q6NZE7-2</a>  |
| Background:             | The function of this protein remains unknown.   |
| Synonyms:               | DKFZp686L20116; MGC131814; RP11-308B5.5   |



[View online »](#)

## Product images:



WB Suggested Anti-Fam122b Antibody  
Titration: 1.0 ug/ml  
Positive Control: Mouse Liver