

Product datasheet for **TA357075**

Otos Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Human, Mouse, Rat Homology: Cow: 86%; Dog: 86%; Guinea Pig: 100%; Human: 100%; Mouse: 100%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	10kDa
Gene Name:	otospiralin
Database Link:	NP_694754 Entrez Gene 260301 Mouse Q8R448
Background:	Otos may be essential for the survival of the neurosensory epithelium of the inner ear.
Synonyms:	OTOSP; otospiralin; OTTHUMP00000164499



[View online »](#)

Product images:



WB Suggested Anti-Otos Antibody
Titration: 1.0 ug/ml
Positive Control: Mouse Pancreas