

Product datasheet for **TA357081**

Nudt16 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Sheep Homology: Cow: 93%; Dog: 100%; Guinea Pig: 100%; Horse: 93%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%; Sheep: 93%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	23kDa
Gene Name:	nudix (nucleoside diphosphate linked moiety X)-type motif 16
Database Link:	NP_083661 Entrez Gene 75686 Mouse Q6P3D0
Background:	The function of this protein remains unknown.
Synonyms:	FLJ31265



[View online »](#)

Product images:



WB Suggested Anti-Nudt16 Antibody
Titration: 1.0 ug/ml
Positive Control: Mouse Small Intestine