

## Product datasheet for **TA357094**

### Fancd2os Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat <b>Homology:</b> Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	19kDa
Gene Name:	FANCD2 opposite strand
Database Link:	<a href="#">NP_001102721</a> <a href="#">Entrez Gene 500294 Rat D3ZY74</a>
Background:	The function of this protein remains unknown.
Synonyms:	RGD1565997



[View online »](#)

## Product images:



Host: Rabbit  
 Target Name: RGD1565997  
 Antibody Dilution: 1.0ug/ml  
 Sample Type: Rat Pancreas