

Product datasheet for **TA357160**

Guk1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Yeast, Zebrafish Homology: Cow: 93%; Dog: 93%; Guinea Pig: 93%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 93%; Rat: 100%; Yeast: 93%; Zebrafish: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	24kDa
Gene Name:	guanylate kinase 1
Database Link:	NP_001013133 Entrez Gene 303179 Rat E9PTV0
Background:	The function of this protein remains unknown.
Synonyms:	FLJ42686; FLJ43710; GMK; OTTHUMP00000037740; OTTHUMP00000037741; OTTHUMP00000037742; OTTHUMP00000037748; OTTHUMP00000037750



[View online »](#)

Product images:



WB Suggested Anti-Guk1 Antibody
Titration: 1.0 ug/ml
Positive Control: Rat Muscle