

Product datasheet for **TA357188**

Gdf10 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat Homology: Cow: 93%; Dog: 93%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 86%; Rat: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	52kDa
Gene Name:	growth differentiation factor 10
Database Link:	NP_852078 Entrez Gene 224650 Mouse Entrez Gene 14560 Mouse P59672
Background:	The function of this protein remains unknown.
Synonyms:	BIP; BMP-3b; BMP3B; GDF-10



[View online »](#)

Product images:



WB Suggested Anti-Gdf10 Antibody
Titration: 1.0 ug/ml
Positive Control: Mouse Brain