

## Product datasheet for **TA357262**

### ANGPTL5 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Dog, Guinea Pig, Horse, Human <b>Homology:</b> Dog: 86%; Guinea Pig: 79%; Horse: 93%; Human: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	44kDa
Gene Name:	angiopoietin like 5
Database Link:	<a href="#">NP_835228</a> <a href="#">Entrez Gene 253935 Human</a> <a href="#">Q86XS5</a>
Background:	The exact functions of ANGPTL5 remain unknown.
Synonyms:	fibrinogen-like
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein



[View online »](#)

## Product images:



WB Suggested Anti-ANGPTL5 Antibody  
Titration: 1.0 ug/ml  
Positive Control: MCF7 Whole Cell