

## Product datasheet for **TA357355**

### Shtn1 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Horse, Human, Mouse, Pig, Rat, Yeast <b>Homology:</b> Horse: 100%; Human: 100%; Mouse: 100%; Pig: 93%; Rat: 100%; Yeast: 85%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	52kDa
Gene Name:	shootin 1
Database Link:	<a href="#">NP_780381</a> <a href="#">Entrez Gene 71653 Mouse</a> <a href="#">Q8K2Q9-2</a>



[View online »](#)

Product images:



WB Suggested Anti-4930506M07Rik Antibody  
Titration: 1.0 ug/ml  
Positive Control: Mouse Heart