

Product datasheet for **TA357425**

Runx1t1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Zebrafish Homology: Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%; Zebrafish: 100%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	65kDa
Gene Name:	runt-related transcription factor 1; translocated to, 1 (cyclin D-related)
Database Link:	NP_001104496 Entrez Gene 12395 Mouse Q8C066
Background:	The function of this protein remains unknown.
Synonyms:	AML1T1; CBFA2T1; CDR; ETO; MGC2796; MTG8; MTG8b; ZMYND2



[View online »](#)

Product images:



WB Suggested Anti-Runx1t1 Antibody
Titration: 1.0 ug/ml
Positive Control: Mouse Brain