

Product datasheet for **TA357474**

Bax Rabbit Polyclonal Antibody


Product data:

| | |
|-------------------------|--|
| Product Type: | Primary Antibodies |
| Applications: | WB |
| Reactivity: | Rat |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | The immunogen is a synthetic peptide directed towards the middle region of Rat Bax |
| Specificity: | Expected reactivity: Cow, Dog, Goat, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Sheep Homology: Cow: 100%; Dog: 100%; Goat: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 93%; Rabbit: 100%; Rat: 100%; Sheep: 100% |
| Formulation: | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i> |
| Concentration: | lot specific |
| Purification: | Affinity Purified |
| Conjugation: | Unconjugated |
| Storage: | For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 21kDa |
| Gene Name: | Bcl2-associated X protein |
| Database Link: | NP_058755 Entrez Gene 24887 Rat G3V8T9 |
| Background: | Bax is a bcl2-related gene and is involved in the regulation of apoptotic cell death. |
| Synonyms: | Bcl2-L-4; BCL2L4 |



[View online »](#)

Product images:



Host: Rabbit
Target Name: Bax
Sample Type: Rat Testis lysates
Antibody Dilution: 1.0ug/ml