

## Product datasheet for **TA357484**

### Atg9a Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Sheep, Zebrafish <b>Homology:</b> Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Rat: 100%; Sheep: 100%; Zebrafish: 86%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	94kDa
Gene Name:	autophagy related 9A
Database Link:	<a href="#">NP_001003917</a> <a href="#">Entrez Gene 245860 Mouse</a> <a href="#">Q3ZAQ4</a>
Background:	The function of this protein remains unknown.
Synonyms:	APG9L1; FLJ22169; MGD3208; OTTHUMP00000206049



[View online »](#)

## Product images:



WB Suggested Anti-Atg9a Antibody  
Titration: 1.0 ug/ml  
Positive Control: Mouse Pancreas