

Product datasheet for TA357509

Product datasneet for 1A557503

AXL Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: WB

Reactivity: Human Host: Rabbit

Clonality: Polyclonal

Immunogen: The immunogen is a synthetic peptide directed towards the C terminal region of human AXL

Specificity: Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Pig, Rat

Homology: Cow: 79%; Dog: 79%; Guinea Pig: 85%; Horse: 93%; Human: 100%; Mouse: 85%;

Pig: 86%; Rat: 85%

Formulation: Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2%

sucrose.

Note that this product is shipped as lyophilized powder to China customers.

Concentration: lot specific

Purification: Affinity Purified
Conjugation: Unconjugated

Storage: For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small

aliquots to prevent freeze-thaw cycles.

Stability: Shelf life: one year from despatch.

Predicted Protein Size: 97kDa

Gene Name: AXL receptor tyrosine kinase

Database Link: NP 001690

Entrez Gene 558 Human

P30530



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



AXL Rabbit Polyclonal Antibody - TA357509


Background: The protein encoded by this gene is a member of the receptor tyrosine kinase subfamily.

Although it is similar to other receptor tyrosine kinases, this protein represents a unique structure of the extracellular region that juxtaposes IgL and FNIII repeats. It transduces signals from the extracellular matrix into the cytoplasm by binding growth factors like vitamin K-dependent protein growth-arrest-specific gene 6. It is involved in the stimulation of cell proliferation and can also mediate cell aggregation by homophilic binding. Alternatively spliced transcript variants encoding different isoforms have been identified.

Synonyms: JTK11; UFO

Protein Families: Druggable Genome, Protein Kinase, Transmembrane

Product images:



WB Suggested Anti-AXL Antibody Titration: 0.2-1

ug/ml

ELISA Titer: 1:312500

Positive Control: Human Placenta