

## Product datasheet for **TA357621**

### OR2M2 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	<b>Expected reactivity:</b> Cow, Dog, Guinea Pig, Horse, Human, Pig, Rat <b>Homology:</b> Cow: 86%; Dog: 86%; Guinea Pig: 79%; Horse: 93%; Human: 100%; Pig: 86%; Rat: 85%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	39kDa
Gene Name:	olfactory receptor family 2 subfamily M member 2
Database Link:	<a href="#">NP_001004688</a> <a href="#">Entrez Gene 391194 Human</a> <a href="#">Q96R28</a>



[View online »](#)

**Background:**

Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms.

**Synonyms:**


OR2M2Q; OST423

**Protein Families:**

Transmembrane

**Protein Pathways:**

Olfactory transduction

**Product images:**

WB Suggested Anti-OR2M2 Antibody  
Titration: 1.0 ug/ml  
Positive Control: THP-1 Whole Cell