

Product datasheet for **TA357649**

Adcyap1r1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity: Cow, Dog, Guinea Pig, Horse, Human, Mouse, Rabbit, Rat, Sheep, Zebrafish Homology: Cow: 100%; Dog: 100%; Guinea Pig: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Rabbit: 93%; Rat: 100%; Sheep: 100%; Zebrafish: 93%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	54kDa
Gene Name:	adenylate cyclase activating polypeptide 1 receptor type 1
Database Link:	NP_598195 Entrez Gene 24167 Rat P32215-2
Background:	Adcyap1r1 regulates neural precursor proliferation; mediates inhibitory signaling for Shh-induced cerebellar granule precursor cell proliferation.
Synonyms:	PAC1; PACAP-R-1; PACAPR; PACAPRI



[View online »](#)

Product images:



WB Suggested Anti-Adcyap1r1 Antibody
Titration: 1.0 ug/ml
Positive Control: Rat Brain