

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA357652

ApInr Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	Expected reactivity : Dog, Horse, Human, Mouse, Rabbit, Rat Homology : Dog: 86%; Horse: 93%; Human: 100%; Mouse: 100%; Rabbit: 86%; Rat: 86%
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. Note that this product is shipped as lyophilized powder to China customers.
Concentration:	lot specific
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	42kDa
Gene Name:	apelin receptor
Database Link:	<u>NP_035914</u> <u>Entrez Gene 23796 Mouse</u> <u>Q9WV08</u>
Background:	Aplnr is a receptor for apelin coupled to G proteins that inhibit adenylate cyclase activity.
Synonyms:	AGTRL1; APJ; APJR; FLJ90771; FLJ96609; HG11; MGC45246



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:



WB Suggested Anti-Aplnr Antibody Titration: 1.0 ug/ml Positive Control: Mouse Pancreas

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US